



Welch Allyn AED 10 Trainer Guidelines 2005 Program

About the Guidelines

On November 28, 2005, American Heart Association released the “2005 Guidelines for Cardiopulmonary Resuscitation and Emergency Cardiovascular Care (ECC)”. Guidelines are updated every five years. The new publication provides several revised emergency care guidelines which include a change to AED protocol. The Welch Allyn AED 10 Trainer upgrade program provides new training protocols, previous Guidelines 2000 protocols and CPR assistance features to support of these new guidelines.

Welch Allyn AED 10 Trainer G5R Program:

AED 10 Trainer *Upgrade* – \$129 (USD)

Program includes:

- Four New G2005 Single-Shock G2005 Training Protocols
- New CPR Coaching Feature With Configurable CPR Compression Rate
 - Welch Allyn’s “Dual Tone Metronome” ensures accurate chest compression and recoil
- Expanded CPR Timer
 - Default 120 seconds (options: 15 through 180)
- Retains 8 Three-Shock G2000 Training Protocols
- New User Manual CD with updated User Manuals in English, French, German, Italian, Spanish, Dutch, and Danish
- Priority processing upon receipt of device
- New set of C-Cell batteries installed

G5R Program Eligibility:

- All Welch Allyn branded AED 10 Trainers are eligible for the upgrade.

To Schedule An Upgrade:

- Contact Welch Allyn’s Technical Service Center at **800-289-2501 – extension 7835** (503.530.7835) to schedule your upgrade.
- An authorization number is required prior to sending the device to our service center, ship only the unit (no accessories such as carrying case, remote, or training pads)
- Coordinator will provide shipping instructions to the Welch Allyn service center

All work done by a certified defibrillation technician. No complicated software to install. Welch Allyn will do all the work for you.